# CHROMagar™ White Opaque Supplement

## Instructions For Use

Available in several languages

**NT-EXT-102**

Version **3.1**

Click below for:

**ENGLISH**

**FRANCAIS**

**ESPAÑOL**

# CHROMagar™ White Opaque Supplement

## MEDIUM PURPOSE

Supplement to be added to the following CHROMagar™ media in order to obtain an opaque background: CHROMagar™ Candida, CHROMagar™ Malassezia, CHROMagar™ O157, CHROMagar™ Orientation, CHROMagar™ Pseudomonas, CHROMagar™ Salmonella Plus, CHROMagar™ Y.enterocolitica, CHROMagar™ mSuperCARBA™, CHROMagar™ LIN-R.

## COMPOSITION

The product is composed of a supplement.

Product	White Opaque Suppl.
Total g/L	1.0 g/L
Aspect	Powder Form
STORAGE	15-30 °C

## PREPARATION (Calculation for 1 L)

### Step 1

Preparation of white opaque background

- Put 1 g of the CHROMagar™ White Opaque supplement reference SU702 in 10 mL sterile water.
- Mix until homogenization. **Appearance: milky suspension without particules/grains.**
- Add the 10 mL to the rehydrated base before the heating step.

### Storage

- Store the prepared media plates in the dark before use.

## PERFORMANCE & LIMITATIONS

This product does not alter the performance of the medium to which is added.

## QUALITY CONTROL

Please perform Quality Control of the final product according to the use of the medium and the local QC regulations and norms.

## WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Inappropriate storage may affect the shelf life of the product.
- Recap the bottles tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.

## DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

## IFU/LABEL INDEX

-  REF Catalogue reference
-  Consult instructions for use
-  Quantity of powder sufficient for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity
-  Protect from light
-  Manufacturer

 Pack Size

Ordering References

5000 mL =

SU702

=

White Opaque Suppl.

25 L =

SU703-25

=

White Opaque Suppl.

CHROMagar™ and Rambach™ are trademarks created by Dr. A. Rambach.  
NT-EXT-102 V3.1 / 07-May-24

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