# Leaflet

LF-EXT-045

Version 5.1

## ● CHROMagar™ C.difficile

Click below:

**EN**

**FR**

**ES**

**CHR**  **Magar™**  
The Chromogenic Media Pioneer

# ● CHROMagar™ C.difficile

**For detection of *Clostridioides difficile***

## Plate Reading

- *C. difficile*  
→ colourless and fluorescent under UV light at 365 nm
- Other bacteria  
→ colourless, not fluorescent or inhibited

## For detection of *Clostridioides difficile*

### Background

*Clostridioides difficile* is the leading cause of nosocomial infectious diarrhea in adults. These infections occur mostly in patients who have both medical care and antibiotic treatment. The symptoms of *C. difficile* infection are fever, abdominal cramps and severe diarrhea leading to death. In the United States, nearly 250,000 people each year develop *C. difficile* infections with at least 14 000 deaths (CDC estimate, 2013). Due to the emergence of highly toxigenic *C. difficile* strains, these infections have become more frequent and more difficult to treat in the last years.

Although PCR has become the leading *C. difficile* detection technique, culture is essential for strain typing and antimicrobial susceptibility testing. CHROMagar™ C.difficile is a fluorogenic culture medium, extremely sensitive and selective, especially designed to simplify and speed up (24 h) the culture of *C. difficile*.

### Intended Use

CHROMagar™ C.difficile is a selective and differential fluorogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with *Clostridioides difficile* to aid in the prevention and control in healthcare settings. The test is performed with stools. Results can be interpreted under UV light after 24 h of anaerobic incubation at 35-37 °C.

CHROMagar™ C.difficile is not intended to diagnose an infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of fluorescent colonies on CHROMagar™ C.difficile does not preclude the presence of *C. difficile*. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

### Medium Performance

1 **RAPID DETECTION COMPARED TO TRADITIONAL MEDIA**  
Big colonies (around 2 mm) of *C. difficile* after only 24 h of incubation in anaerobic atmosphere, contrary to traditional media requiring 48 h.

2 **HIGH SENSITIVITY AND SPECIFICITY**

	Analytical data *	Clinical data**	
		CHROMagar™ C. difficile	Reference medium (Taurocholate-CCFA)
Sensitivity	-	95.4 %	70 %
Specificity	100 %	88.8 %	97 %

\* Data obtained after a 24 h incubation at 35 °C under anaerobic atmosphere in the study «Comparison of CHROMagar C.difficile and taurocholate-CCFA media for isolation of toxigenic *Clostridium difficile* from stools». CCFA = Cefoxitin-Cycloserine-Fructose Agar. Roux et al. Poster ASM 2014.

\*\* Data obtained after a 24 h incubation at 35 °C under anaerobic atmosphere with 594 stool samples, being positive 174, in the study «Comparison of CHROMagar C.difficile and taurocholate-CCFA media for isolation of toxigenic *Clostridium difficile* from stools». Roux et al. Poster ASM 2014.

### Medium Description

<b>Powder Base</b>	Total ..... 54.7 g/L Agar ..... 15.0 Peptone & Yeast extract ..... 25.0 Salts ..... 9.0 Growth factors ..... 4.0 Chromogenic mix ..... 1.7 Storage at 15/30 °C - pH: 7.8 ± 0.2 Shelf Life ..... > 18 months
<b>Supplement</b> (Included in the pack)	Powder form ..... 325 mg/L Storage at 2/8 °C ..... Shelf Life ..... > 12 months
Usual Samples	Stools
Procedure	Direct streaking. Incubation at 35-37 °C for 24 h. Anaerobic condition.
Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com	

### Order References

Please use these product references when contacting your local distributor:

5000 mL pack ..... CD122  
(included in this reference: powder base CD122(B) + supplement CD122(S))

Manufacturer: CHROMagar, 29 avenue George Sand,  
93210 La Plaine Saint-Denis - France  
Email: CHROMagar@CHROMagar.com  
Website: www.CHROMagar.com  
Find your nearest distributor on  
www.chromagar.com/en/find-your-distributor/

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