CHROMagar™ Candida


Click below:

EN

FR

ES

● CHROMagar™ Candida

**For isolation and differentiation
of major clinical-significant *Candida* species**




Plate Reading

• *Candida albicans*
→ green

• *Candida tropicalis*
→ metallic blue

• *Candida krusei*
→ pink, fuzzy

For isolation and differentiation of major clinical-significant *Candida* species

Background

Perhaps as a result of progress in medical therapy with a related decrease in patient resistance and increase in the number of susceptible hosts, the incidence of bloodstream *Candida* infections rose through the 1990s but has stabilised in recent years. More commonly, *Candida* species are involved in superficial oropharyngeal and urogenital infections particularly for immuno-depressed people such as the elderly and AIDS victims. "Early diagnosis is essential for early effective management of the patients." (WHO Guidelines on Standard Operating Procedures for Laboratory Diagnosis of HIV-Opportunistic Infections)

Although *C. albicans* remains the major species involved, other types such as *C. tropicalis*, *C. krusei* or *C. glabrata* have increased proportionately as new antifungal agents have worked very effectively against *C. albicans*. This shows the importance of an accurate detection for a proper antifungal therapy choice.

Intended Use

CHROMagar™ Candida is a selective chromogenic culture medium intended for use in the qualitative direct detection, differentiation and presumptive identification of *Candida* species. The test is performed with swabs from skin, throat, ears and vaginal specimens as well as sputum, urine and stools samples, in parallel to cultures on Sabouraud agar, to aid in the Candidiasis diagnosis. Results can be interpreted after 20-48 h of aerobic incubation at 35-37 °C.

Further microbiological testing or epidemiological typing are needed.

A lack of growth or the absence of colonies on CHROMagar™ Candida does not preclude the presence of *Candida*.

CHROMagar™ Candida is not intended to diagnose infection nor to guide nor monitor treatment for infections.

Medium Performance

- 1 **HIGH CONTRAST**
Intense colony colouration helping to differentiate species and mixed cultures.
- 2 **FAST**
20/48 h incubation.
- 3 **PANORAMIC VIEW**
Panoramic view of the *Candida* population mix with the ability to detect major species but also the presence of minor populations.
- 4 **POWERFUL TOOL**

	Analytical data *	Clinical data **
	Sensitivity & Specificity	CHROMagar™ Candida Sensitivity & Specificity
<i>C. albicans</i>	100 %	96.6 % / 97.9 %
<i>C. tropicalis</i>	97.9 % / 98.8 %	100 % / 98.5 %
<i>C. krusei</i>	100 %	-
<i>C. glabrata</i>	97.7 % / 98.4 %	100 %

* Data obtained after a 48 h incubation at 30 °C in aerobic conditions in the study "A comparison of methods for yeast identification including CHROMagar Candida, Vitek system YBC and a traditional bio chemical method". Huang *et al.*, 2001. *Chinese Med. J.*

** Data obtained after a 48 h incubation at 37 °C in aerobic conditions with 127 clinical specimens from patients in the study "Evaluation of chromogenic media and seminested PCR in the identification of *Candida* species". Daef *et al.*, 2014. *Brazilian J. Microbiol.*

- 5 **LESS WORKLOAD**
With traditional media such as Sabouraud agar, identification of mixed cultures is difficult, whereas CHROMagar™ Candida, just by the colour of the colony, will instantly differentiate various *Candida* species and help locating their colonies in mixed cultures.

Medium Description

Powder Base	Total
Total	47.7 g/L
Agar	15.0
Peptone	10.2
Chromogenic mix	22.0
Chloramphenicol	0.5
Storage at 15/30 °C - pH: 6.1 ± 0.2.	
Shelf Life	> 18 months

Usual Samples	Swabs from skin, throat, ears and vaginal specimens, sputum, urine and stools.
Procedure	Direct Streaking, Incubation 20-48 h, 35-37 °C. Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Order References

Please use these product references when contacting your local distributor:

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