
CHROMagar™ StrepB


Click below:

EN

FR

ES

● CHROMagar™ StrepB

**For the isolation and differentiation of
Group B *Streptococcus* (*S. agalactiae*)**

For the isolation and differentiation of Group B *Streptococcus* (*S. agalactiae*)

Background

GBS in pregnant women: The Group B *Streptococci* (GBS), also known under the name of *Streptococcus agalactiae*, are the cause of numerous infections in adults but mostly an important cause of serious neonatal infections, occurring in the three first weeks of life. Studies indicate that approximately 12-27 % of pregnant women are colonised by GBS. (WHO, Infectious diseases, Group B *Streptococcus*). Detecting vaginal (and in some countries also rectal) colonisation by GBS in pregnant women is the most effective strategy to prevent infection transmission during baby delivery.

Worldwide, official guidelines recommend prenatal screening of GBS in the last month of pregnancy. In GBS-carrying women, this screening allows determining the need of intrapartum antibiotic prophylaxis, which has been proved effective in preventing the infections occurring in the first week of life (known as early-onset GBS infections).

Intended Use

CHROMagar™ StrepB is a selective chromogenic culture medium that is intended to aid in the qualitative determination of Group B *Streptococcus* (GBS) colonization in pregnant women. This medium supports the growth of hemolytic and non-hemolytic GBS strains. The test is performed by direct, or after enrichment, plating of urine, vaginal and/or vagino/rectal swabs obtained from pregnant women. CHROMagar™ StrepB results can be interpreted after 18-24 h aerobic incubation at 35-37 °C.

CHROMagar™ StrepB is not intended to diagnose infection nor to guide nor monitor treatment for infections. Further identification, susceptibility testing and epidemiological typing is needed on suspect colonies.

Medium Performance

Compared to CNA Blood Agar and Granada:

- 1 **EASY INTERPRETATION**
Easier reading thanks to an intense mauve colony colouration.
- 2 **HIGH SENSITIVITY**
Detection of GBS, including non-haemolytic strains, with sensitivity close to 100 %.
- 3 **HIGH SPECIFICITY**
Differentiation of GBS from other bacteria by selective inhibition or by counter-colouration.
- 4 **FAST**
Results in 18-24 h.
- 5 **SIMPLICITY**
Incubation in aerobic conditions. No need of CO₂. Confirmation by latex agglutination can be performed directly from the colony.

	Analytical data *	Clinical data**	
		CHROMagar™ StrepB	Reference medium (CNA)
Sensitivity	100 %	94 %	92 %
Specificity	80 %	100 %	100 %

* In-house data obtained after 24 h incubation at 37 °C in aerobic conditions.

** Data obtained after 24 h incubation at 35 °C in aerobic conditions with 242 vaginal and rectal swabs in the study «Evaluation of four chromogenic media for the isolation of Group B *Streptococcus* from vaginal specimens in pregnant women». Salem & Anderson, 2015. *Pathology*.



Plate Reading

- Group B *Streptococcus*
→ mauve
- Other Microorganism
→ blue, colourless or inhibited

Medium Description

Powder Base	Total 44.7 g/L Agar 15.0 Peptones and yeast extract 20.0 Salts 7.5 Chromogenic mix 2.2 Storage at 15/30 °C - pH: 7.3 +/- 0.2 Shelf Life > 18 months
+	
2 Supplements (included in the pack)	1 st : Liquid form..... 8 mL/L 2 nd : Powder..... 0.25 g/L Storage at 15/30 °C Storage at 2/8 °C Shelf Life > 18 months Shelf Life > 18 months

Usual Samples	Vaginal samples, rectal swabs, urine
Procedure	Direct Streaking. Possible prior enrichment with our LIM RambaQUICK™ StrepB (ref. LB082). Incubation 18-24 h at 35-37 °C. Aerobic conditions.
Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com	


Order References

Please use these product references when contacting your local distributor:

- 5000 mL pack SB282
(included in this reference: powder base SB282(B) + supplements SB282(S1) + SB282(S2))
25 L pack SB283-25
(included in this reference: powder base SB283-25(B) + supplements SB283-25(S1) + SB283-25(S2))

Manufacturer: CHROMagar, 29 avenue George Sand, 93210 La Plaine Saint-Denis - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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