● CHROMagar™ ESBL

Click below:

EN

FR

ES

● CHROMagar™ ESBL

**For overnight detection of Gram-negative bacteria
producing Extended Spectrum Beta-Lactamase**

For overnight detection of Gram-negative bacteria producing Extended Spectrum Beta-Lactamase

Background

ESBL (Extended Spectrum Beta-Lactamases) are enzymes that mediate resistance to penicillins, extended-spectrum third generation cephalosporins (C3G) and monobactams. ESBL-producing Enterobacteriaceae started to appear in the 80's, first emerging in *K. pneumoniae* and recognized through significant hospital outbreaks worldwide. In the early 2000's, ESBL have settled in the commensal species *Escherichia coli* and spread into the human community as an asymptomatic colonizer. But patient can become infected from this reservoir, and the early detection of ESBL-producing bacteria carriers is important to minimize the spread of colonization in their environment.

Intended Use

CHROMagar™ ESBL is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with Extended-spectrum-beta-lactamase-resistant Enterobacteria (ESBL) to aid in the prevention and control of ESBL in healthcare settings. The test is performed with rectal swab and stools from patients to screen for ESBL colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C. The medium can also be used as an early warning indicator for diagnostic tests of infections to signal the possible presence of multi drug-resistant bacteria. This use does not replace the institution's protocols.

CHROMagar™ ESBL is not intended to diagnose infection caused by ESBL producers nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ ESBL does not preclude the presence of ESBL-producing bacteria. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

Medium Performance

1 FAST RESULTS
Detection after overnight incubation

2 SPECIES DIFFERENTIATION
thanks to the chromogenic performances of supplemented CHROMagar™ Orientation. Indeed, the product is composed of a powder base CHROMagar™ Orientation and a supplement to enhance ESBL-producing bacteria.

3 HIGH SENSITIVITY AND SPECIFICITY

	Analytical data *		Clinical data **	
	CHROMagar™ ESBL	CHROMagar™ ESBL	CHROMagar™ ESBL	Competitor's medium
Sensitivity	100 %	98.3 %	97.5 %	
Specificity	97 %	72.3 %	72.9 %	

* Data obtained after a 18-24 h incubation at 37 °C in aerobic conditions in the study « Rapid detection of ESBL-producing Enterobacteria on chromogenic medium: Colorex Orientation/ESBL ». Laudat et al. Poster SFM 2010.

** Data obtained by testing 2 337 rectal screening samples, on plates incubated at 37 °C for 18 h, in the study « Comparison of five media for detection of Extended-Spectrum beta-Lactamase by use of the wrap instrument for automated specimen processing ». Grohs et al., 2013. *J.Clin. Microbiol.*

4 TIME AND WORKLOAD SAVINGS
Direct culture from specimen. There is no need of a selective pre-enrichment.

5 FLEXIBILITY
CHROMagar™ ESBL is supplied with a more than 18 months shelf-life . This allows flexibility of use, whether in an epidemic situation with many patients to screen, or whether for random surveillance of cultures.


Plate Reading

- *E. coli* ESBL
→ dark pink to reddish
- *Klebsiella*, *Enterobacter*, *Citrobacter* ESBL
→ metallic blue (+/- red halo)
- *Proteus* ESBL
→ brown halo
- *Pseudomonas* ESBL
→ translucent cream to blue
- *Acinetobacter* ESBL
→ cream, opaque

Medium Description

Powder Base (CHROMagar™ Orientation)	Total 33 g/L Agar 15.0 Peptone and yeast extract 17.0 Chromogenic mix 1.0 Storage at 15/30 °C - pH: 7.0 +/-0.2 Shelf Life > 18 months
+	
CHROMagar™ ESBL supplement (included in the pack)	Selective mix (Powder form) 0.57 g/L Storage at 2/8 °C Shelf Life > 18 months

Usual Samples	Rectal swabs and stools
Procedure	Direct Streaking. Incubation 18-24h at 35-37 °C. Aerobic conditions

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Order References

Please use these product references when contacting your local distributor:


- 5000 mL pack ESRT2 (included in this reference: powder base RT412 + supplement ES372)
- 25 L pack ESRT3-25 (included in this reference: powder base RT413-25 + supplement ES373-25)

Manufacturer: CHROMagar, 29 avenue George Sand,
93210 La Plaine Saint-Denis - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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