
## CHROMagar™ MH Orientation


Click below:

[EN](#)

[FR](#)

[ES](#)

# CHROMagar™ MH Orientation

**Chromogenic Mueller Hinton medium**

## Chromogenic Mueller Hinton medium

### Background

CHROMagar™ MH Orientation is a chromogenic Mueller Hinton agar. This product was developed with the aim to speed up the time to obtain results by combining the chromogenic differentiation of the species in the sample, and the antibiotic susceptibility test. This tool allows a same-day results in both fields, the identification and susceptibility, directly in the primary culture, contrary to the concomitant traditional protocol which takes 48 hours.

An evaluation was carried by AB Biodisk in 2008, comparing the performances of our CHROMagar™ MH Orientation with the traditional Mueller-Hinton formulation<sup>1</sup>. This study, carried with over 20 classical antimicrobial agents, concluded that the use of our medium with E-test strips can provide same day results.

CHROMagar™ MH Orientation has been used principally in two applications, but it certainly may find usefulness in many other cases.

The first application is on samples from ICU patients with Ventilated Associated Pneumonia (VAP), which is one of the most frequent ICU nosocomial infections with high morbidity and mortality. Shortening the period in which empirical therapy can be optimized results in better outcomes. As stated by the Spanish team, leader in this protocol<sup>2</sup>, plating respiratory samples directly onto CHROMagar™ MH Orientation is "a rapid procedure for antimicrobial susceptibility testing (...) crucial for modifying therapeutic regimens".

The second application is the common Urine Tract Infections (UTI) for which physicians usually prescribe an empiric therapy. Urine samples with typical infection counts of 10<sup>3</sup>-10<sup>5</sup> cfu/mL can be spread directly on this plate. Pluri-microbial samples as well as monomicrobial-contaminated samples can be easily investigated by this method.

In this case also as in the VAP, the same-day identification and susceptibility results can confirm the empirical therapy, or help in the decision to change it.

#### References:

<sup>1</sup> Engelhardt A. *et al*, Rapid E-test® MIC Testing Using CHROMagar™ And Mueller Hinton Agar For Gram Positive And Gram Negative Aerobes - ECCMID 2008

<sup>2</sup> E. Cercenado *et al*, Evaluation of direct E-test on lower respiratory tract samples using a chromogenic agar medium: a rapid procedure for antimicrobial susceptibility testing - ECCMID 2009

### Medium Performance

#### 1 2 IN 1

It allows chromogenic identification and antibiotic susceptibility in the same plate.

#### 2 FAST RESULTS

For many microorganisms/antibiotics, the results can be read just after overnight incubation.


#### 3 EASY PLATE READING

Thanks to the chromogenic reactions, diameters and/or inhibitions zones are clearly defined.

### Medium Description


<b>Powder Base</b>	Total ..... 39.4 g/L Agar ..... 17.0 Peptones & growth factors ..... 21.0 Chromogenic mix ..... 1.4 Storage at 15/30 °C - pH: 7.1 ± 0.2 Shelf Life ..... > 18 months
<b>+</b>	
<b>Supplement</b> (included in the pack)	Growth and regulator factors ..... 4 mL/L Storage at 15/30 °C Aspect: Liquid Form Shelf Life ..... > 18 months

Usual Samples	urine, respiratory excretions
Procedure	Direct spreading. Incubation 18-24 h at 35-37 °C. Aerobic conditions.
For detailed preparation procedure, please refer to our IFU.	


### Antibiogram Plate Reading

- *E. coli*  
→ dark pink to reddish
- *Klebsiella*, *Enterobacter*, *Citrobacter*  
→ metallic blue


#### UTI Test Analysis: on CHROMagar™ MH Orientation

--> *Klebsiella* + *Pseudomonas*

#### UTI Test Analysis: on CHROMagar™ MH Orientation

--> *E. coli* + *Enterococci*  
(contamination)

#### E-test !: --> mix of Gram(-)

**E-test !:**  
*E. faecalis* after 8 h incubation on basic MH agar

**E-test !:**  
*E. faecalis* after 8 h incubation on CHROMagar™ MH Orientation agar

## Order References


Please use these product references when contacting your local distributor:

- Pack de 5000 mL ..... MH482  
(included in this reference: powder base MH482(B) + supplement MH482(S))
- Pack de 25 L ..... MH483-25  
(included in this reference: powder base MH483-25(B) + supplement 5 x MH482(S))

Manufacturer: CHROMagar, 29 avenue George Sand, 93210 La Plaine Saint-Denis - France  
Email: CHROMagar@CHROMagar.com  
Website: www.CHROMagar.com

Find your nearest distributor on  
www.CHROMagar.com/contact

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<sup>2</sup> E. Cercenado et al, Evaluation of direct E-test on lower respiratory tract samples using a chromogenic agar medium: a rapid procedure for antimicrobial susceptibility testing - ECCMID 2009

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