## CHROMagar™ C3G<sup>R</sup>

Click below:

# CHROMagar™ C3G<sup>R</sup>

**For rapid detection (18-24h) of beta-lactamase producing  
Enterobacteria achieving third generation cephalosporin  
resistance**

## For overnight detection of beta-lactamase producing Enterobacteria achieving third generation cephalosporin resistance

### Background

β-Lactamase production (ESBL, AmpC, carbapenemase) is the most common mechanism of β-lactam drug resistance in Gram negative bacteria. Each of the involved mechanisms can be plasmid-encoded and transfer horizontally. Moreover, plasmids continuously gather other antibiotic resistance determinants, ultimately closing to pan-resistance. Their spread leads to increased colonization pressure in hospital and community as well. It is crucial to ensure that proper surveillance is in place to help implement appropriate guidelines and policies for their spread control. Rapid detection of bacteria producing these enzymes also allows for de-escalation to more targeted therapy, to conserve carbapenem antibiotics for more serious infections.

### Intended Use

CHROMagar™ C3G<sup>R</sup> is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with 3<sup>rd</sup> generation cephalosporin-resistant Enterobacteria (C3G<sup>R</sup>-E) to aid in the prevention and control of C3G<sup>R</sup>-E in healthcare settings. The test is performed with rectal swab and stools from patients to screen for C3G<sup>R</sup>-E colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

The medium can also be used as an early warning indicator for diagnostic tests of infections to signal the possible presence of multi drug-resistant bacteria. This use does not replace the institution's protocols. CHROMagar™ C3G<sup>R</sup> is not intended to diagnose C3G<sup>R</sup>-E infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ C3G<sup>R</sup> does not preclude the presence of C3G<sup>R</sup>-E. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.


### Medium Performance

- 1 **FAST RESULTS**  
Detection after overnight incubation.
- 2 **SPECIES DIFFERENTIATION**  
Thanks to the chromogenic performances of supplemented CHROMagar™ Orientation. Indeed, the product is composed of a powder base CHROMagar™ Orientation and a supplement to select β-Lactamase producing bacteria.
- 3 **HIGH SENSITIVITY**  
Unique medium not inhibiting plasmid-mediated AmpC-producing bacteria.
- 4 **TIME AND WORKLOAD SAVINGS**  
Direct culture from specimen. There is no need of a selective pre-enrichment.



	Analytical data *	Clinical data**	
		CHROMagar™ C3G <sup>R</sup>	Reference medium (ESBL Isolation Agar Non chromogenic media)
Sensitivity	100 %	77 %	67 %
Specificity	-	-	-

\* Data obtained after 24 h incubation at 37 °C in aerobic conditions in the study «Evaluation of Four Media for the Isolation of Resistant Gram Negative Organisms from Surveillance Rectal Swabs». Joshi-Caesar *et al.* Poster ASM 2012.

\*\* Data obtained after 24 h incubation at 37 °C in aerobic conditions with 96 rectal samples in the study «Evaluation of Four Media for the Isolation of Resistant Gram Negative Organisms from Surveillance Rectal Swabs». Joshi-Caesar *et al.* Poster ASM 2012.

### Plate Reading

-  • *E. coli*  
→ dark pink to reddish
-  • *Klebsiella, Enterobacter, Citrobacter*  
→ metallic blue (+/- red halo)
-  • *Proteus*  
→ brown halo
-  • *Pseudomonas*  
→ translucent cream to blue
-  • *Acinetobacter*  
→ cream, opaque

### Medium Description

Powder Base (CHROMagar™ Orientation)	Total .....	33 g/L
	Agar .....	15.0
	Peptone and yeast extract .....	17.0
	Chromogenic mix .....	1.0
	Storage at 15/30 °C - pH: 7.0 +/- 0.2	
	Shelf Life .....	> 18 months
+ CHROMagar™ C3G <sup>R</sup> Supplement (included in the pack)	Selective mix (Powder form) .....	0.37 g/L
	Storage at 2/8 °C	
	Shelf Life .....	> 18 months

Usual Samples	stools and rectal specimens.
Procedure	Direct Streaking. Incubation 18-24 h at 35-37 °C. Aerobic conditions

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document) available on [www.CHROMagar.com](http://www.CHROMagar.com)

### Order References

Please use these product references when contacting your local distributor:

- 5000 mL pack ..... CGRT2  
(included in this reference: powder base RT412 + supplement CG632)
- 25 L pack ..... CGRT3-25  
(included in this reference: powder base RT413-25 + supplement CG633-25)

Manufacturer: CHROMagar, 29 avenue George Sand,  
93210 La Plaine Saint-Denis - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)  
Find your nearest distributor on  
[www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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