# Leaflet

LF-EXT-057  
Version 11.1

## CHROMagar™ mSuperCARBA™


Click below:


CHR  Magar™  
The Chromogenic Media Pioneer




# CHROMagar™ mSuperCARBA™

**For detection and isolation of Carbapenem-Resistant  
*Enterobacteriaceae* (CRE)**

## Plate Reading

-  • CPE *E. coli*  
→ dark pink to reddish
-  • CPE *Klebsiella, Enterobacter, Citrobacter*  
→ metallic blue
-  • CPO\* *Pseudomonas*  
→ translucent cream to blue
-  • CPO\* *Acinetobacter*  
→ cream, opaque

\*CPE : Carbapenemase-producing enterobacteriaceae

\*CPO : Carbapenemase-producing organism

## For detection and isolation of Carbapenem-resistant *Enterobacteriaceae* (CRE)

### Background

CDC: «Carbapenem-resistant *Enterobacteriaceae* (CRE) are usually resistant to all  $\beta$ -lactam agents as well as most other classes of antimicrobial agents. The treatment options for patients infected with CRE are very limited. Healthcare-associated outbreaks of CRE have been reported. Patients colonized with CRE are thought to be a source of transmission in the healthcare setting. Identifying patients who are colonized with CRE and placing these patients in isolation precautions may be an important step in preventing transmission».

CHROMagar™ launched in 2007 the first chromogenic medium for the detection of carbapenem-resistant bacteria, particularly targeting KPC-enzymes. Since then, many other Carbapenemase-producing organisms (CPO) have been spreading around the world and therefore there was a need to address today the difficult detection of low level CPO.

Alain Rambach and Patrice Nordmann have joined their efforts to develop a highly sensitive chromogenic medium, CHROMagar™ mSuperCARBA™, the new generation of chromogenic media that reaches unprecedented performances: detection of a large variety of CPO : KPC, NDM, VIM, IMP, OXA...with an impressive limit of detection (10 CFU/mL) even for weakly expressed CPO like OXA-48, while maintaining a high level of selectivity.

### Intended Use

CHROMagar™ mSuperCARBA™ is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with carbapenem-resistant *Enterobacteria* (CRE), including OXA-48 producers, to aid in the prevention and control of CRE in healthcare settings. The test is performed with rectal swab and stools from patients to screen for CRE colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

CHROMagar™ mSuperCARBA™ is not intended to diagnose CRE infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ mSuperCARBA™ does not preclude the presence of CRE. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

### Medium Performance

- 1 HIGHLY SENSITIVE**  
Most CPO (including OXA-48 and OXA-48 like) detected after overnight incubation.
- 2 IMPRESSIVE LIMIT OF DETECTION**  
10 CFU/mL
- 3 HIGHLY SELECTIVE AND SPECIFIC**  
Not only beta-lactam susceptible bacteria are inhibited, but also most of the ESBL and AmpC hyperproducers are inhibited, rendering the tool very specific for carbapenemase-producing bacteria.

### Medium Description

<b>Powder Base</b>	Total ..... 42.5 g/L Agar ..... 15.0 Peptones ..... 20.0 Salt ..... 5.0 Chromogenic and selective mix ..... 0.8 Growth factors ..... 1.7 Storage at 15/30 °C - pH: 7.2 +/- 0.2 Shelf Life ..... > 18 months
<b>2 Supplements (included in the pack)</b>	1 <sup>st</sup> : Liquid form.....2 mL/L Storage at 15/30 °C Shelf Life ..... > 18 months
	2 <sup>nd</sup> : Powder..... 0.25 g/L Storage at 2/8 °C Shelf Life ..... > 18 months

<b>Usual Samples</b>	Rectal swabs and stools.
<b>Procedure</b>	Direct Streaking, Incubation 18-24 h at 35-37 °C Aerobic conditions.

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document) available on [www.CHROMagar.com](http://www.CHROMagar.com)

	Analytical data *	Clinical Data **
	CHROMagar™ mSuperCARBA™	
Sensitivity	100 %	100 %
Specificity	71 %	100 %

\* Data obtained after 24 h incubation at 37 °C in aerobic conditions in the study «Amélioration de la détection des Entérobactéries Productrices de Carbapénémase (EPC)». Dos Santos *et al.* RICA 2017.

\*\* Data obtained after 24 h incubation at 35 °C with 211 rectal swabs from the study «CHROMagar™ mSuperCARBA: performance in carbapenem-resistant *Enterobacteriaceae* isolates characterized at molecular level and routine surveillance rectal swabs specimens». García-Fernández *et al.* 2017. *Diagn. Microbiol. Infect. Dis.*

### Order References

Please use these product references when contacting your local distributor:

5000 mL pack ..... SC172 (included in this reference: powder base SC172) + suppl. SC172(S1) + SC172(S2))

25 L pack ..... SC173-25

(included in this reference: powder base SC173-25) + suppl. SC173-25(S1) + SC173-25(S2))

Manufacturer: CHROMagar, 29 avenue George Sand, 93210 La Plaine Saint-Denis - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)

Find your nearest distributor on [www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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