
Leaflet

LF-EXT-060

Version 5.1

CHROMagar™ COL-APSE

Click below:


EN

FR

ES

CHR  Magar™
The Chromogenic Media Pioneer

● CHROMagar™ COL-APSE

**For detection of Colistin resistant
Gram-negative bacteria**


Plate Reading

- COL-R *E. coli*
→ Dark pink to reddish
- COL-R *Klebsiella*, *Enterobacter*, *Citrobacter*
→ Metallic blue
- COL-R *Pseudomonas*
→ Translucent cream to green
- COL-R *Acinetobacter*
→ Cream, opaque

For detection of Colistin resistant Gram-negative bacteria

Background

Polymyxin E (colistin) and B are increasingly used as antimicrobials in the treatment of multi-drug resistant bacterial infections. Polymyxin resistance, although intrinsic in Gram-positive and some Gram-negative species (*Proteus*, *Morganella*, *Serratia*), is now a problem in a number of other pathogens (*Acinetobacter baumannii*, *Pseudomonas aeruginosa*, *Escherichia coli*, *Salmonella enterica*, *Klebsiella pneumoniae*).¹

Resistance arises due to mutations/insertions in genes involved in LPS biosynthesis (lpx, pmrA/B, mgrB, phoP/Q) and/or the acquisition of phosphoethanolamine transferases (PEtN). Of great concern is the recently described plasmid-encoded PEtN, MCR-1, now found worldwide in a range of animal, human and environmental bacterial isolates.¹

CHROMagar™ COL-APSE is a sensitive and specific medium for the growth of Colistin resistant bacterial pathogens with a lower limit of detection of 10 CFU/mL. This new chromogenic medium may be useful as a primary isolation medium in the surveillance and recovery of Colistin resistant bacteria from complex human, veterinary and environmental samples especially those with plasmid mediated MCR-1 or novel mechanisms of polymyxin resistance.¹

1- Novel Chromogenic Culture Media (CHROMagar™ COL-APSE) for the Isolation and Differentiation of Colistin Resistant Gram-negative Pathogens. Muhd Haziq F Abdul Momin, David W Wareham. ECCMID 2017.

Intended Use

CHROMagar™ COL-APSE is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with colistin-resistant Gram-negative bacteria (COL-R) to aid in the prevention and control of COL-R in healthcare settings. The test is performed with rectal swab, perineal swabs and stools samples from patients to screen for COL-R colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

The medium can also be used as an early warning indicator for diagnostic tests of infections to signal the possible presence of COL-R bacteria. This use does not replace the institution's protocols. CHROMagar™ COL-APSE is not intended to diagnose COL-R infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ COL-APSE does not preclude the presence of COL-R. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

CHROMagar™ COL-APSE can also be used in the detection of COL-R in the analyses of food products for human consumption, animal feed, in livestock and in environmental samples.

Medium Performance

1 COLOURFUL DIFFERENTIATION OF COLONIES WITH ACQUIRED COLISTIN RESISTANCE

Distinguishing between *E. coli*, Coliforms, *Pseudomonas* and *Acinetobacter*.

2 IMPRESSIVE LIMIT OF DETECTION

10 CFU/mL

	Analytical data *	Clinical data **
	CHROMagar™ COL-APSE	
Sensitivity	100 %	-
Specificity	81 %	-
NPV	-	100 %
PPV	-	83 %

* Data obtained after a 24 h incubation at 37 °C in aerobic conditions in the study «CHROMagar COL-APSE: a selective bacterial culture medium for the isolation and differentiation of colistin-resistant Gram-negative pathogens». Abdul Momin et al., 2017. *J. Med. Microbiol.*

** Data obtained after a 24 h incubation at 37 °C in aerobic conditions with 98 human stool samples in the study «Wide dissemination of colistin-resistant *Escherichia coli* with the mobile resistance gene mcr in healthy residents in Vietnam». Yamamoto et al., 2018. *J. Antimicrob. Chemother.*

Medium Description

Powder Base	Total	42.5 g/L
	Agar	15.0
	Peptones	20.0
	Salt	5.0
	Chromogenic and selective mix	0.8
	Growth factors	1.7
	Storage at 15/30 °C - pH : 7.1 +/-0.2	
	Shelf Life	> 18 months
+		
Supplement (included in the pack)	Liquid form.....	2 mL/L
	Storage at 15/30 °C	
	Shelf Life	> 18 months

Usual Samples	rectal swab, perineal swab, stools, food products, animal feed, livestock and environmental samples.
Procedure	Direct Streaking, Incubation 18-24 h at 35-37 °C Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Order References

Please use these product references when contacting your local distributor:

5000 mL pack CO262
(included in this reference: powder base CO262(B) + supplement CO262(S))


Manufacturer: CHROMagar, 29 avenue George Sand,
93210 La Plaine Saint-Denis - France

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