CHROMagar™ StrepA


Click below:

EN

FR

ES

● CHROMagar™ StrepA

**For screening of Group A Streptococci
in throat samples**

For screening of Group A Streptococci in throat samples

Background

Streptococcus pyogenes (GAS) is responsible for 5 % to 15 % consultations of pediatricians and other physicians for sore throat, with approximately 15 million consultations per year in the United States (Shulman *et al.*, 2012). However, the most common cause of acute pharyngitis remains viral and does not require antibiotic treatment.

GAS affects only 20-30 % of pharyngitis in children and 5-15 % in adults, but antibiotic treatment is critical to avoid many invasive and immunopathological complications. The decision for this treatment must often be made in a consulting room or in an emergency department. Several clinical scores have been developed, but it is often impossible to distinguish the streptococcal etiology by clinical aspects alone.

For this reason, diagnostic algorithms have been developed involving tests that can be performed during the consultation or that confirm the prescriber's decision as quickly as possible: Antigenic tests, molecular tests and culture on blood agar, which remains the historical Gold Standard. Only the culture allows to perform an antibiotic susceptibility test and the epidemiological follow-up of the strains.

Intended Use

CHROMagar™ StrepA is a selective chromogenic culture medium intended for use in the qualitative direct detection, differentiation and presumptive identification of *Streptococcus pyogenes* -group A *Streptococcus* (GAS) to aid in the diagnosis of GAS throat infections. The test is performed with throat swabs. Results can be interpreted after 18-24 h of incubation at 35-37 °C under CO₂ atmosphere.

Concomitant cultures are necessary to recover organisms for further microbiological testing or epidemiological typing. A lack of growth or the absence of colonies on CHROMagar™ StrepA does not preclude the presence of GAS. CHROMagar™ StrepA is not intended to diagnose infection nor to guide nor monitor treatment for infections.

Medium Performance

- 1 **FIRST CHROMOGENIC MEDIUM**
Designed for the detection of *S. pyogenes*.
- 2 **EASY INTERPRETATION**
Easier reading thanks to an intense orange to red colony colouration.
- 3 **HIGH SENSITIVITY AND SPECIFICITY**


	Analytical data *	Clinical data**
Sensitivity	100 %	96.7 %
Specificity	97.5 %	100 %

* Data obtained after a 24 h incubation at 35-37 °C in aerobic conditions in the study "Digital detection and the use of artificial intelligence to detect Group A *Streptococcus* using a chromogenic agar". Keyak *et al.*, Poster ASM 2019.

** Data obtained after a 20 h incubation at 35 °C in a CO₂ atmosphere with 159 throat samples being positive 120, in the study "Validation and implementation of Colorex™ CHROMagar™ Strep A agar on WASP™/WASPLab™ for screening for *Streptococcus pyogenes* using the ESwab™". Gaskin *et al.*, Poster ASM 2019.

Plate Reading

- Group A *Streptococcus*
→ orange to red
- Other oral Streptococci
→ colourless or blue

Medium Description

Powder Base	Total 54.3 g/L Agar 15.0 Peptones and yeast extract 27.0 Salts 9.0 Growth factors 3.0 Chromogenic and selective mix 0.3 Storage at 15/30 °C - pH: 7.2 +/- 0.2 Shelf Life > 12 months
+	
2 Supplements (included in the pack)	1 st : Liquid form 2 mL/L 2 nd : Powder 0.1 g/L Storage at 15/30 °C Storage at 2/8 °C Shelf Life > 18 months Shelf Life > 18 months

Usual Samples	Throat swabs
Procedure	Direct Streaking. Incubation 18-24 h at 35-37 °C. Under CO ₂ .

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Order References

Please use these product references when contacting your local distributor:

- 5000 mL pack SP372
(included in this reference: powder base SP372(B) + supplements SP372(S1) and SP372(S2))
- 25 L pack SP373-25
(included in this reference: powder base SP373-25(B) + supplements SP373-25(S1) and SP373-25(S2))
- Bulk on request

Manufacturer: CHROMagar, 29 avenue George Sand, 93210 La Plaine Saint-Denis - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com

Find your nearest distributor on www.CHROMagar.com/contact

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