Leaflet

LF-EXT-069

Version 2.1

● CHROMagar™ Serratia

Click below:

EN

FR

ES

CHR  Magar™
The Chromogenic Media Pioneer

● CHROMagar™ Serratia

For detection of *Serratia marcescens*



Plate Reading

• *S. marcescens*
→ Green-blue to metallic blue

Medium Description

| | |
|---|--|
| Powder Base | Total 42.5 g/L Agar 15.0 Peptones 20.0 Salt 5.0 Growth factors 1.7 Chromogenic and selective mix 0.8 Storage at 15/30 °C - pH: 7.1 +/- 0.2 Shelf Life > 18 months |
| + | |
| Supplement (included in the pack) | Liquid form 2 mL/L Storage at 15/30 °C Shelf Life > 18 months |
| Usual Samples | Rectal, throat, wound swabs, stools and surface samples. |
| Procedure | Direct Streaking. Incubation 18-24 h at 35-37 °C Aerobic conditions. |
| Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com | |

For detection of *Serratia marcescens*

Background

Serratia species are implicated in nosocomial infections. In several countries, *Serratia marcescens* is frequently associated with epidemics in intensive care units and in particular in neonatal and pediatric units. Surveillance of nosocomial infections requires effective recovery of clinical isolates from faeces, wound exudates and respiratory samples to prevent problems of cross infection and potentially fatal infections. In this context, CHROMagar™ Serratia, a culture medium perfectly suited to the search for *S. marcescens* in faeces.

In addition, *S. marcescens* is able to survive days to months on surfaces, distilled water and hand soap, making it an important pathogen in nosocomial infections and sporadic epidemics.

Intended Use

CHROMagar™ Serratia is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of colonization with *Serratia marcescens* to aid in the prevention and control of *Serratia marcescens* in intensive care unit (ICU), notably neo-natal. The test is performed with rectal swabs, throat swabs, wound swabs and stools from patients to screen for *Serratia* colonization. It can also be used in hygiene monitoring in the clinical environment with surface sampling. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

CHROMagar™ Serratia is not intended to diagnose *Serratia* infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ Serratia does not preclude the presence of *Serratia*. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

Medium Performance

- RELIABLE**
First commercially available chromogenic medium for *Serratia marcescens*.
- VERY HIGH SENSITIVITY AND SELECTIVITY**

| | Analytical data * | Clinical data** |
|-------------|---------------------|-----------------|
| | CHROMagar™ Serratia | |
| Sensitivity | 100 % | 97 % |
| Specificity | 97 % | 100 % |

* Data obtained after 20 h incubation at 37 °C in aerobic conditions in the study «Validation of Colorex™ (CHROMagar™) Serratia agar on WASP™/WASPLab™ in screening for *Serratia marcescens* in neonatal intensive care units using the ESwab™». Gaskin *et al.* Poster ECCMID 2020.

** Data obtained after 24 h incubation at 37 °C in aerobic conditions with 96 epidemiological samples (46 rectal swabs and 50 from environmental surfaces, sink and drains) in the study «Evaluation of CHROMagar™ Serratia agar, a new chromogenic medium for the detection and isolation of *Serratia marcescens*». Pérez Viso *et al.*, 2021. *Eur. J. Clin. Microbiol. Infect. Dis.*

- HIGH SPECIFICITY AND LOWER WORKLOAD**
The conventional media used for the detection of *S. marcescens* like MacConkey Agar has very poor specificity, creating an abundance of false positives. Because all lactose negative bacteria have a similar aspect than *S. marcescens*. On the contrary, CHROMagar™ Serratia detect all *Serratia* (those pigmented or not) and allows technicians to focus on the real contaminated samples.
- INTENSE GREEN-BLUE COLOURATION FOR EASY READING**
compared to conventional media (MacConkey, Blood Agar).
- FAST RESULTS**
particularly useful in case of sudden outbreak of *Serratia marcescens*.

Order References


Please use these product references when contacting your local distributor:

5000 mL pack SM302 (included in this reference: powder base SM302(B) + suppl. SM302(S))

Manufacturer: CHROMagar, 29 avenue George Sand, 93210 La Plaine Saint-Denis - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com

Find your nearest distributor on www.CHROMagar.com/contact

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