
# Leaflet

LF-EXT-073

Version 3.1

# CHROMagar™ Candida Plus

Click below:


**EN**

**FR**

**ES**

**CHR**  **Magar™**  
The Chromogenic Media Pioneer

# CHROMagar™ Candida Plus

**For detection and differentiation of major clinical *Candida* species, including *C. auris***

## For detection and differentiation of major clinical *Candida* species, including *C. auris*

### Background

*Candida* are yeast species involved in various infections called candidiasis, which can affect damaged skin, respiratory tract, digestive and urogenital systems. These candidiasis can be severe with significant morbidity for nosocomial infections or in immunocompromised patients. Although *C. albicans* is still the main species involved, the use of antifungal agents has given rise to other species such as *C. tropicalis*, *C. krusei* and *C. glabrata*.

In 2016, the World Health Organization added *C. auris* to this list, with a prevalence of over 90 % resistant to fluconazole. In addition, some strains are multidrug resistant to amphotericin B, voriconazole, and/or echinocandins.

CHROMagar™ Candida Plus is the first chromogenic isolation medium to detect and differentiate *C. auris* in addition to other major clinical *Candida* species such as *C. albicans*, *C. tropicalis*, *C. glabrata* or *C. krusei*.

### Intended Use

CHROMagar™ Candida Plus is a selective chromogenic culture medium intended for use in the qualitative direct detection, differentiation and presumptive identification of *Candida* species including *C. auris*. The test is performed with swabs from skin, throat, ears and vaginal specimens as well as sputum, urine and stools samples, in parallel to cultures on Sabouraud agar, to aid in the Candidiasis diagnosis. Results can be interpreted after 24-48 h of aerobic incubation at 30-37 °C.

Further microbiological testing or epidemiological typing are needed. A lack of growth or the absence of colonies on CHROMagar™ Candida Plus does not preclude the presence of *Candida*.

CHROMagar™ Candida Plus is not intended to diagnose infection nor to guide nor monitor treatment for infections.

### Medium Performance

#### 1 HIGH SPECIFICITY

Differentiation of the most common *Candida* species with very high specificity (analytical data):

<i>C. albicans</i>	100 % *
<i>C. tropicalis</i>	100 % *
<i>C. krusei</i>	100 % *

#### 2 UNIQUE MEDIUM TO DIFFERENTIATE *C. auris* FROM OTHER *Candida* SPECIES

Owe to its high specificity, it can also be used also as a screening tool in case of outbreaks for patient samples and surfaces suspected of being contaminated with *C. auris*

For <i>C. auris</i> :	Specificity	100 % **
	Sensitivity	100 % **

\* Data obtained after a 48 h incubation at 37 °C in aerobic conditions in the study « Evaluation of a novel chromogenic medium for *Candida* spp. identification and comparison with CHROMagar™ Candida for the detection of *Candida auris* in surveillance samples. Mulet Bayona et al., 2020. *Diag. Microbiol. Inf. Dis.*

\*\* Data obtained after a 24-48 h incubation at 37 °C in aerobic conditions with 364 patients surveillance samples and 212 environmental samples in the study «Novel chromogenic medium CHROMagar™ Candida Plus for detection of *Candida auris* and other *Candida* species from surveillance and environmental samples: A multicenter study. Mulet Bayona et al., 2022. *J. of Fungi.*

#### 3 EASY IDENTIFICATION



Identification by MALDI-TOF can be carried directly from a colony. No need of subculture.

### Medium Description

Powder Base	Total .....	50.9 g/L
	Agar .....	15.0
	Peptones .....	11.0
	Chromogenic and selective mix .....	24.9
	Storage at 15/30 °C - pH: 6.1 ± 0.2.	
	Shelf Life .....	> 9 months

Usual Samples	Swabs from skin, throat, ears, vagina as well as sputum, urine and stools samples.
Procedure	Direct Streaking. Incubation 24-48 h, 30-37 °C. Aerobic conditions.

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document) available on [www.CHROMagar.com](http://www.CHROMagar.com)

### Plate Reading

- *C. auris*  
→ Light blue with blue halo  
Blue from the back side
- *C. albicans*  
→ Green-blue
- *C. tropicalis*  
→ Metallic blue with pink halo
- *C. glabrata*  
→ Mauve
- *C. krusei*  
→ Pink and fuzzy

### Order References

Please use these product references when contacting your local distributor:

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